

ZeroG-Passive is an easy-to-use static body-weight support system that is compatible with the ZeroG Gait and Balance System track. ZeroG-Passive provides patients the ability to safely practice overground walking and balance activities. The patient's harness connects to the Passive trolley that glides smoothly along the ceiling track.

The patient is protected against falls as the Passive trolley maintains its orientation above the patient. Because the trolley weighs less than 5 pounds and uses precision bearings, it is nearly imperceptible to patients. With ZeroG-Passive, patients can begin walking early and intensively after neurological and orthopedic injuries, factors related to enhanced outcomes.

